









Reusable Spo2 Sensor

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2151-060A	03-0101-379	Mindray	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	
M2251-060A	03-0101-380	MEC-1000/1200	One-Piece Adult Soft Tip Spo2 Sensor		
M2151-060E	03-0101-381	PM7000/8000/9000	One-Piece Adult Ear Clip Spo2 Sensor		
M2151-060P	03-0101-382	0	One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-060P	03-0101-383	PM600,509B	One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-060I	03-0101-384	300cm Long	One-Piece Infant Soft Tip Spo2 Sensor		
M2351-060W	03-0101-385		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-060Y	03-0101-387		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-066A	03-0101-388	Mindray	One-Piece Adult Finger Clip Spo2 Sensor	 7Pin Assembled Connector	
M2251-066A	03-0101-389	PM6800,iPM	One-Piece Adult Soft Tip Spo2 Sensor		
M2151-066E	03-0101-390	Beneview T5/T6/T8	One-Piece Adult Ear Clip Spo2 Sensor		
M2151-066P	03-0101-391	300cm Long	One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-066P	03-0101-392		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-066I	03-0101-393		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-066W	03-0101-394		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-066Y	03-0101-396		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-034A	03-0101-541	Mindray	One-Piece Adult Finger Clip Spo2 Sensor	 12Pin Molded Connector	
M2251-034A	03-0101-542	PM5000/6000	One-Piece Adult Soft Tip Spo2 Sensor		
M2151-034E		300cm Long	One-Piece Adult Ear Clip Spo2 Sensor		
M2151-034P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-034P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-034I			One-Piece Infant Soft Tip Spo2 Sensor		
M2351-034W			One-Piece Neonatal Wrap Spo2 Sensor		
M2351-034Y			One-Piece Y-Style Multi-use Spo2 Sensor		



SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2151-059AM	03-0101-397	Mindray-Masimo PM9000 Express PM6000/7000 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	       
M2251-059AM	03-0101-398		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-059EM	03-0101-399		One-Piece Adult Ear Clip Spo2 Sensor		
M2151-059PM	03-0101-400		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-059PM	03-0101-401		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-059IM	03-0101-402		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-059WM	03-0101-403		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-059YM	03-0101-405		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-060AM	03-0101-406	Mindray-Masimo PM9000 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	
M2251-060AM	03-0101-407		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-060EM	03-0101-408		One-Piece Adult Ear Clip Spo2 Sensor		
M2151-060PM	03-0101-409		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-060PM	03-0101-410		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-060IM	03-0101-411		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-060WM	03-0101-412		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-060YM	03-0101-414		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-034AM		Mindray-Masimo PM6000 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 12Pin Molded Connector	
M2251-034AM	03-0101-543		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-034EM			One-Piece Adult Ear Clip Spo2 Sensor		
M2151-034PM			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-034PM			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-034IM			One-Piece Infant Soft Tip Spo2 Sensor		
M2351-034WM			One-Piece Neonatal Wrap Spo2 Sensor		
M2351-034YM			One-Piece Y-Style Multi-use Spo2 Sensor		




SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2151-076A	03-0101-415	Mindray-Masimo BeneView T5/T8 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 7Pin Assembled Connector	       
M2251-076A	03-0101-416		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-076E	03-0101-417		One-Piece Adult Ear Clip Spo2 Sensor		
M2151-076P	03-0101-418		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-076P	03-0101-419		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-076I	03-0101-420		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-076W	03-0101-421		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-076Y	03-0101-423		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-087AC	03-0101-433	Mindray-Masimo With Chip iMEC 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 8Pin Assembled Connector	    
M2251-087AC	03-0101-434		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-087EC	03-0101-435		One-Piece Adult Ear Clip Spo2 Sensor		
M2151-087PC	03-0101-436		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-087PC	03-0101-437		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-087IC	03-0101-438		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-087WC	03-0101-439		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-087YC	03-0101-441		One-Piece Y-Style Multi-use Spo2 Sensor		
M2151-067A	03-0101-424	Mindray- Oximax BeneView T5/T8 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 8Pin Assembled Connector	  
M2251-067A	03-0101-425		One-Piece Adult Soft Tip Spo2 Sensor		
M2151-067E	03-0101-426		One-Piece Adult Ear Clip Spo2 Sensor		
M2151-067P	03-0101-427		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2251-067P	03-0101-428		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2251-067I	03-0101-429		One-Piece Infant Soft Tip Spo2 Sensor		
M2351-067W	03-0101-430		One-Piece Neonatal Wrap Spo2 Sensor		
M2351-067Y	03-0101-432		One-Piece Y-Style Multi-use Spo2 Sensor		


SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2152-063A	03-0101-442	Goldway UT4000A/B/F Choice... 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 5Pin Assembled Connector	       
M2252-063A	03-0101-443		One-Piece Adult Soft Tip Spo2 Sensor		
M2152-063E	03-0101-444		One-Piece Adult Ear Clip Spo2 Sensor		
M2152-063P	03-0101-445		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2252-063P	03-0101-446		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2252-063I	03-0101-447		One-Piece Infant Soft Tip Spo2 Sensor		
M2352-063W	03-0101-448		One-Piece Neonatal Wrap Spo2 Sensor		
M2352-063Y	03-0101-450		One-Piece Y-Style Multi-use Spo2 Sensor		
M1153-001AD		Comen- Digital DB9 90cm Long	Adult Finger Clip Spo2 Sensor	 DB9 Molded Connector	
M1253-001AD			Adult Soft Tip Spo2 Sensor		
M1153-001ED			Adult Ear Clip Spo2 Sensor		
M1153-001PD			Pediatric Finger Clip Spo2 Sensor		
M1253-001PD			Pediatric Soft Tip Spo2 Sensor		
M1253-001ID			Infant Soft Tip Spo2 Sensor		
M1353-001WD			Neonatal Wrap Spo2 Sensor		
M1353-001YD			Y-Style Multi-use Spo2 Sensor		
M2153-092AC	03-0101-564	Comen-Masimo With Chip 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 12Pin Molded Connector	
M2253-092AC			One-Piece Adult Soft Tip Spo2 Sensor		
M2153-092EC			One-Piece Adult Ear Clip Spo2 Sensor		
M2153-092PC			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2253-092PC			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2253-092IC			One-Piece Infant Soft Tip Spo2 Sensor		
M2353-092WC			One-Piece Neonatal Wrap Spo2 Sensor		
M2353-092YC			One-Piece Y-Style Multi-use Spo2 Sensor		




SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2153-061A	03-0101-562	Comen Zondan 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	       
M2253-061A			One-Piece Adult Soft Tip Spo2 Sensor		
M2153-061E			One-Piece Adult Ear Clip Spo2 Sensor		
M2153-061P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2253-061P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2253-061I			One-Piece Infant Soft Tip Spo2 Sensor		
M2353-061W			One-Piece Neonatal Wrap Spo2 Sensor		
M2353-061Y			One-Piece Y-Style Multi-use Spo2 Sensor		
M2153-061AD	03-0101-505	Comen-Digital Start-8000A/B/C Start-8000D/H Start-5000B/C 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	
M2253-061AD	03-0101-506		One-Piece Adult Soft Tip Spo2 Sensor		
M2153-061ED	03-0101-507		One-Piece Adult Ear Clip Spo2 Sensor		
M2153-061PD	03-0101-508		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2253-061PD	03-0101-509		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2253-061ID	03-0101-510		One-Piece Infant Soft Tip Spo2 Sensor		
M2353-061WD	03-0101-511		One-Piece Neonatal Wrap Spo2 Sensor		
M2353-061YD	03-0101-513		One-Piece Y-Style Multi-use Spo2 Sensor		
M2154-061A		Guoteng GT3000/5000 GT6000/7000 GT9000 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	
M2254-061A			One-Piece Adult Soft Tip Spo2 Sensor		
M2154-061E			One-Piece Adult Ear Clip Spo2 Sensor		
M2154-061P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2254-061P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2254-061I			One-Piece Infant Soft Tip Spo2 Sensor		
M2354-061W			One-Piece Neonatal Wrap Spo2 Sensor		
M2354-061Y			One-Piece Y-Style Multi-use Spo2 Sensor		


SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M1156-001AD	03-0101-451	Biolight-Digital 90cm Long	Adult Finger Clip Spo2 Sensor	 DB9 Molded Connector	       
M1256-001AD	03-0101-452		Adult Soft Tip Spo2 Sensor		
M1156-001ED	03-0101-453		Adult Ear Clip Spo2 Sensor		
M1156-001PD	03-0101-454		Pediatric Finger Clip Spo2 Sensor		
M1256-001PD	03-0101-455		Pediatric Soft Tip Spo2 Sensor		
M1256-001ID	03-0101-456		Infant Soft Tip Spo2 Sensor		
M1356-001WD	03-0101-457		Neonatal Wrap Spo2 Sensor		
M1356-001YD	03-0101-459		Y-Style Multi-use Spo2 Sensor		
M2156-063A	03-0101-460	Biolight M6,M12,M69-Old 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 5Pin Assembled Connector	
M2256-063A	03-0101-461		One-Piece Adult Soft Tip Spo2 Sensor		
M2156-063E	03-0101-462		One-Piece Adult Ear Clip Spo2 Sensor		
M2156-063P	03-0101-463		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2256-063P	03-0101-464		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2256-063I	03-0101-465		One-Piece Infant Soft Tip Spo2 Sensor		
M2356-063W	03-0101-466		One-Piece Neonatal Wrap Spo2 Sensor		
M2356-063Y	03-0101-468		One-Piece Y-Style Multi-use Spo2 Sensor		
M2156-063AD	03-0101-469	Biolight-Digital M8000,M8000A M9000,M9000A BDT-352A,M69 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 5Pin Assembled Connector	
M2256-063AD	03-0101-470		One-Piece Adult Soft Tip Spo2 Sensor		
M2156-063ED	03-0101-471		One-Piece Adult Ear Clip Spo2 Sensor		
M2156-063PD	03-0101-472		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2256-063PD	03-0101-473		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2256-063ID	03-0101-474		One-Piece Infant Soft Tip Spo2 Sensor		
M2356-063WD	03-0101-475		One-Piece Neonatal Wrap Spo2 Sensor		
M2356-063YD	03-0101-477		One-Piece Y-Style Multi-use Spo2 Sensor		



SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2156-077AD	03-0101-478	Biolight-Digital M900,M7000 M9500,A8 Modular 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 12Pin Assembled Connector	       
M2256-077AD	03-0101-479		One-Piece Adult Soft Tip Spo2 Sensor		
M2156-077ED	03-0101-480		One-Piece Adult Ear Clip Spo2 Sensor		
M2156-077PD	03-0101-481		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2256-077PD	03-0101-482		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2256-077ID	03-0101-483		One-Piece Infant Soft Tip Spo2 Sensor		
M2356-077WD	03-0101-484		One-Piece Neonatal Wrap Spo2 Sensor		
M2356-077YD	03-0101-486		One-Piece Y-Style Multi-use Spo2 Sensor		
M2156-083AD	03-0101-487	Biolight-Digital A8 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 9Pin Assembled Connector	
M2256-083AD	03-0101-488		One-Piece Adult Soft Tip Spo2 Sensor		
M2156-083ED	03-0101-489		One-Piece Adult Ear Clip Spo2 Sensor		
M2156-083PD	03-0101-490		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2256-083PD	03-0101-491		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2256-083ID	03-0101-492		One-Piece Infant Soft Tip Spo2 Sensor		
M2356-083WD	03-0101-493		One-Piece Neonatal Wrap Spo2 Sensor		
M2356-083YD	03-0101-495		One-Piece Y-Style Multi-use Spo2 Sensor		
M2156-084AD	03-0101-496	Biolight-Digital V6 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 7Pin Assembled Connector	
M2256-084AD	03-0101-497		One-Piece Adult Soft Tip Spo2 Sensor		
M2156-084ED	03-0101-498		One-Piece Adult Ear Clip Spo2 Sensor		
M2156-084PD	03-0101-499		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2256-084PD	03-0101-500		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2256-084ID	03-0101-501		One-Piece Infant Soft Tip Spo2 Sensor		
M2356-084WD	03-0101-502		One-Piece Neonatal Wrap Spo2 Sensor		
M2356-084YD	03-0101-504		One-Piece Y-Style Multi-use Spo2 Sensor		



SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2157-060AD		Tomorrow- Digital 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	       
M2257-060AD			One-Piece Adult Soft Tip Spo2 Sensor		
M2157-060ED			One-Piece Adult Ear Clip Spo2 Sensor		
M2157-060PD			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2257-060PD			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2257-060ID			One-Piece Infant Soft Tip Spo2 Sensor		
M2357-060WD			One-Piece Neonatal Wrap Spo2 Sensor		
M2357-060YD			One-Piece Y-Style Multi-use Spo2 Sensor		
M2157-060A			Tomorrow 300cm Long		
M2257-060A		One-Piece Adult Soft Tip Spo2 Sensor			
M2157-060E		One-Piece Adult Ear Clip Spo2 Sensor			
M2157-060P		One-Piece Pediatric Finger Clip Spo2 Sensor			
M2257-060P		One-Piece Pediatric Soft Tip Spo2 Sensor			
M2257-060I		One-Piece Infant Soft Tip Spo2 Sensor			
M2357-060W		One-Piece Neonatal Wrap Spo2 Sensor			
M2357-060Y		One-Piece Y-Style Multi-use Spo2 Sensor			
M2157-061AD		Tomorrow- Digital 300cm Long		One-Piece Adult Finger Clip Spo2 Sensor	
M2257-061AD			One-Piece Adult Soft Tip Spo2 Sensor		
M2157-061ED			One-Piece Adult Ear Clip Spo2 Sensor		
M2157-061PD			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2257-061PD			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2257-061ID			One-Piece Infant Soft Tip Spo2 Sensor		
M2357-061WD			One-Piece Neonatal Wrap Spo2 Sensor		
M2357-061YD			One-Piece Y-Style Multi-use Spo2 Sensor		

SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2159-060AD	03-0101-544	General-Digital G2A,G3A,G3C G3D,G3F,G3G G3H,G3L,G6A... 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	       
M2259-060AD			One-Piece Adult Soft Tip Spo2 Sensor		
M2159-060ED			One-Piece Adult Ear Clip Spo2 Sensor		
M2159-060PD			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2259-060PD			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2259-060ID			One-Piece Infant Soft Tip Spo2 Sensor		
M2359-060WD			One-Piece Neonatal Wrap Spo2 Sensor		
M2359-060YD			One-Piece Y-Style Multi-use Spo2 Sensor		
M1163-001AD	03-0101-514		Create-Digital DB9 300cm Long		
M1263-001AD	03-0101-515	Adult Soft Tip Spo2 Sensor			
M1163-001ED	03-0101-516	Adult Ear Clip Spo2 Sensor			
M1163-001PD	03-0101-517	Pediatric Finger Clip Spo2 Sensor			
M1263-001PD	03-0101-518	Pediatric Soft Tip Spo2 Sensor			
M1263-001ID	03-0101-519	Infant Soft Tip Spo2 Sensor			
M1363-001WD	03-0101-520	Neonatal Wrap Spo2 Sensor			
M1363-001YD	03-0101-522	Y-Style Multi-use Spo2 Sensor			
M2163-062AD	03-0101-523	Create-Digital UP-7000,UP-9000 PC-900A/S,PC-5000A/B NewTech:VS V Neu Vision500 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 5Pin Assembled Connector	   
M2263-062AD	03-0101-524		One-Piece Adult Soft Tip Spo2 Sensor		
M2163-062ED	03-0101-525		One-Piece Adult Ear Clip Spo2 Sensor		
M2163-062PD	03-0101-526		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2263-062PD	03-0101-527		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2263-062ID	03-0101-528		One-Piece Infant Soft Tip Spo2 Sensor		
M2363-062WD	03-0101-529		One-Piece Neonatal Wrap Spo2 Sensor		
M2363-062YD	03-0101-531		One-Piece Y-Style Multi-use Spo2 Sensor		

SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2164-061A	03-0101-532	Choice 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	       
M2264-061A	03-0101-533		One-Piece Adult Soft Tip Spo2 Sensor		
M2164-061E	03-0101-534		One-Piece Adult Ear Clip Spo2 Sensor		
M2164-061P	03-0101-535		One-Piece Pediatric Finger Clip Spo2 Sensor		
M2264-061P	03-0101-536		One-Piece Pediatric Soft Tip Spo2 Sensor		
M2264-061I	03-0101-537		One-Piece Infant Soft Tip Spo2 Sensor		
M2364-061W	03-0101-538		One-Piece Neonatal Wrap Spo2 Sensor		
M2364-061Y	03-0101-540		One-Piece Y-Style Multi-use Spo2 Sensor		
M2168-063A		Yonke 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 5Pin Assembled Connector	
M2268-063A			One-Piece Adult Soft Tip Spo2 Sensor		
M2168-063E			One-Piece Adult Ear Clip Spo2 Sensor		
M2168-063P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2268-063P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2268-063I			One-Piece Infant Soft Tip Spo2 Sensor		
M2368-063W			One-Piece Neonatal Wrap Spo2 Sensor		
M2368-063Y			One-Piece Y-Style Multi-use Spo2 Sensor		
M2169-061A	03-0101-549	Adecon 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	 6Pin Assembled Connector	
M2269-061A			One-Piece Adult Soft Tip Spo2 Sensor		
M2169-061E			One-Piece Adult Ear Clip Spo2 Sensor		
M2169-061P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2269-061P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2169-061I			One-Piece Infant Soft Tip Spo2 Sensor		
M2369-061W			One-Piece Neonatal Wrap Spo2 Sensor		
M2369-061W			One-Piece Y-Style Multi-use Spo2 Sensor		

SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M1101-300A		2-Lead LED 660nm/ 905nm For:Nellcor, BCI Philips..... 90cm Long	Adult Finger Clip Spo2 Sensor	Without Connector	       
M1201-300A			Adult Soft Tip Spo2 Sensor		
M1101-300E			Adult Ear Clip Spo2 Sensor		
M1101-300P			Pediatric Finger Clip Spo2 Sensor		
M1201-300P			Pediatric Soft Tip Spo2 Sensor		
M1201-300I			Infant Soft Tip Spo2 Sensor		
M1301-300W			Neonatal Wrap Spo2 Sensor		
M1301-300Y			Y-Style Multi-use Spo2 Sensor		
M2101-300A		2-Lead LED 660nm/ 905nm For:Nellcor, BCI Philips..... 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	Without Connector	
M2201-300A			One-Piece Adult Soft Tip Spo2 Sensor		
M2101-300E			One-Piece Adult Ear Clip Spo2 Sensor		
M2101-300P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2201-300P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2201-300I			One-Piece Infant Soft Tip Spo2 Sensor		
M2301-300W			One-Piece Neonatal Wrap Spo2 Sensor		
M2301-300Y			One-Piece Y-Style Multi-use Spo2 Sensor		
M2105-300A		Digital-235/237	One-Piece Adult Finger Clip Spo2 Sensor	Without Connector	
M2205-300A		2-Lead LED 660nm/ 905nm For:Dolphin Biolight..... 300cm Long	One-Piece Adult Soft Tip Spo2 Sensor		
M2105-300E			One-Piece Adult Ear Clip Spo2 Sensor		
M2105-300P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2205-300P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2205-300I			One-Piece Infant Soft Tip Spo2 Sensor		
M2305-300W			One-Piece Neonatal Wrap Spo2 Sensor		
M2305-300Y			One-Piece Y-Style Multi-use Spo2 Sensor		

SPO2

Model	Part Number	Original Brand	Description	Sensor Connector	Sensor Photo
M2103-300A		3-Lead LED 660nm/ 940nm For:CSI, Ohmeda Nihon Konden..... 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	Without Connector	       
M2203-300A			One-Piece Adult Soft Tip Spo2 Sensor		
M2103-300E			One-Piece Adult Ear Clip Spo2 Sensor		
M2103-300P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2203-300P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2203-300I			One-Piece Infant Soft Tip Spo2 Sensor		
M2303-300W			One-Piece Neonatal Wrap Spo2 Sensor		
M2303-300Y			One-Piece Y-Style Multi-use Spo2 Sensor		
M2106-300A		3-Lead LED 660/ 880/ 940nm For:Dixtal/Brazil 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	Without Connector	
M2206-300A			One-Piece Adult Soft Tip Spo2 Sensor		
M2106-300E			One-Piece Adult Ear Clip Spo2 Sensor		
M2106-300P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2206-300P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2206-300I			One-Piece Infant Soft Tip Spo2 Sensor		
M2306-300W			One-Piece Neonatal Wrap Spo2 Sensor		
M2306-300Y			One-Piece Y-Style Multi-use Spo2 Sensor		
M2107-300A		4-Lead LED 660/ 880/ 940nm For:GE-OxyTip 300cm Long	One-Piece Adult Finger Clip Spo2 Sensor	Without Connector	
M2207-300A			One-Piece Adult Soft Tip Spo2 Sensor		
M2107-300E			One-Piece Adult Ear Clip Spo2 Sensor		
M2107-300P			One-Piece Pediatric Finger Clip Spo2 Sensor		
M2207-300P			One-Piece Pediatric Soft Tip Spo2 Sensor		
M2207-300I			One-Piece Infant Soft Tip Spo2 Sensor		
M2307-300W			One-Piece Neonatal Wrap Spo2 Sensor		
M2307-300Y			One-Piece Y-Style Multi-use Spo2 Sensor		

SPO2