



### Reusable NIBP Cuff (PU Cuff with TPU Bladder)

P/N	Size	Colour	Model	Description	Weight&Package	Photo
03-0201-033	Infant	Orange	M5111/Single	Cuff Width:6.5cm	25g	 <p style="text-align: center;">Tube Material:PVC Tube Length:20cm</p>
03-0201-035			M5211/Double	Cuff Length:21cm	35g	
03-0201-037	Pediatric	Green	M5112/Single	Limb Circumference:9-14cm	1pcs/Bag	
03-0201-039			M5212/Double	Cuff Width:9.5cm	40g	
03-0201-041	Small Adult	Light Blue	M5113/Single	Cuff Length:30cm	50g	
03-0201-043			M5213/Double	Limb Circumference:14-21cm	1pcs/Bag	
03-0201-045	Adult	Brown	M5114/Single	Cuff Width:11.5cm	55g	
03-0201-047			M5214/Double	Cuff Length:39.5cm	65g	
03-0201-049	Adult Extra	Brown	M5114L/Single	Limb Circumference:21-27cm	1pcs/Bag	
03-0201-051			M5214L/Double	Cuff Width:14.5cm	90g	
03-0201-053	Adult Extra Large	Brown	M5114XL/Single	Cuff Length:50cm	100g	
03-0201-055			M5214XL/Double	Limb Circumference:27-35cm	100g	
03-0201-057	Large Adult	Deep Red	M5115/Single	Limb Circumference:27-35/44cm	1pcs/Bag	
03-0201-059			M5215/Double	Cuff Width:15.5cm	100g	
03-0201-061	Thigh	Gray	M5116/Single	Cuff Length:56cm	110g	
03-0201-063			M5216/Double	Limb Circumference:27-35/53cm	115g	
				Cuff Width:17.5cm	115g	
				Cuff Length:55cm	125g	
				Limb Circumference:35-44cm	1pcs/Bag	
				Cuff Width:21.5cm	170g	
				Cuff Length:74.5cm	180g	
				Limb Circumference:44-53cm	1pcs/Bag	

NIBP

## Reusable NIBP Cuff


(Nylon Laminated TPU Cuff With Bladderless)

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-097	Neonate	White	M5301/Single	Cuff Width:5.7cm	25g	 <p>Tube Material:PVC Tube Length:20cm</p>
03-0201-099			M5401/Double	Cuff Length:21.5cm Limb Circumference:9-14.5cm	35g 1pcs/Bag	
03-0201-101	Infant	White	M5302/Single	Cuff Width:8.2cm	35g	
03-0201-103			M5402/Double	Cuff Length:30.5cm Limb Circumference:13-21.5cm	45g 1pcs/Bag	
03-0201-105	Pediatric	White	M5303/Single	Cuff Width:10.5cm	50g	
03-0201-107			M5403/Double	Cuff Length:41.5cm Limb Circumference:21-27cm	60g 1pcs/Bag	
03-0201-109	Adult	White	M5304/Single	Cuff Width:13.5cm	70g	
03-0201-111			M5404/Double	Cuff Length:51.5cm Limb Circumference:26-35.5cm	80g 1pcs/Bag	
03-0201-113	Adult Extra	White	M5304L/Single	Cuff Width:13.5cm	90g	
03-0201-115			M5404L/Double	Cuff Length:66cm Limb Circumference:26-35.5/45cm	100g 1pcs/Bag	
03-0201-117	Adult Extra Large	White	M5304XL/Single	Cuff Width:14cm	110g	
03-0201-119			M5404XL/Double	Cuff Length:720cm Limb Circumference:26-35.5/55cm	120g 1pcs/Bag	
03-0201-121	Large Adult	White	M5305/Single	Cuff Width:16.6cm	95g	
03-0201-123			M5405/Double	Cuff Length:66.5cm Limb Circumference:35-45cm	105g 1pcs/Bag	
03-0201-125	Thigh	White	M5306/Single	Cuff Width:20.6cm	125g	
03-0201-127			M5406/Double	Cuff Length:79cm Limb Circumference:44-55cm	135g 1pcs/Bag	

NIBP

## Reusable NIBP Cuff


(Polyester Cuff With TPU Bladder)

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-001	Neonate	White	M5121/Single	Cuff Width:6cm	20g	 <p style="text-align: center;">Tube Material:PVC Tube Length:20cm</p>
03-0201-003			M5221/Double	Cuff Length:21cm	30g	
03-0201-005	Infant	White	M5122/Single	Limb Circumference:6-11cm	1pcs/Bag	
03-0201-007			M5222/Double	Cuff Width:8.8cm	30g	
03-0201-009	Pediatric	White	M5123/Single	Cuff Length:29cm	40g	
03-0201-011			M5223/Double	Limb Circumference:10-19cm	1pcs/Bag	
03-0201-013	Adult	White	M5124/Single	Cuff Width:11.5cm	50g	
03-0201-015			M5224/Double	Cuff Length:39cm	60g	
03-0201-017	Adult Extra	White	M5124L/Single	Limb Circumference:18-26cm	1pcs/Bag	
03-0201-019			M5224L/Double	Cuff Width:14.5cm	70g	
03-0201-021	Adult Extra Large	White	M5124XL/Single	Cuff Length:50cm	80g	
03-0201-023			M5224XL/Double	Limb Circumference:25-35/44cm	80g	
03-0201-025	Large Adult	White	M5125/Single	Cuff Length:56cm	90g	
03-0201-027			M5225/Double	Limb Circumference:25-35/53cm	110g	
03-0201-029	Thigh	White	M5126/Single	Cuff Length:69cm	120g	
03-0201-031			M5226/Double	Limb Circumference:34-47cm	160g	
				Limb Circumference:46-66cm	170g	

NIBP

## Reusable NIBP Cuff

(Nylon Cuff With TPU Bladder)

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-065	Neonate	White	M5101/Single	Cuff Width:6cm	20g	 <p>Tube Material:PVC Tube Length:20cm</p>
03-0201-067			M5201/Double	Cuff Length:21cm	30g	
03-0201-069	Infant	White	M5102/Single	Limb Circumference:6-11cm	1pcs/Bag	
03-0201-071			M5202/Double	Cuff Width:8.8cm	30g	
03-0201-073	Pediatric	White	M5103/Single	Cuff Length:29cm	40g	
03-0201-075			M5203/Double	Limb Circumference:10-19cm	1pcs/Bag	
03-0201-077	Adult	White	M5104/Single	Cuff Width:11.5cm	50g	
03-0201-079			M5204/Double	Cuff Length:39cm	60g	
03-0201-081	Adult Extra	White	M5104L/Single	Limb Circumference:18-26cm	1pcs/Bag	
03-0201-083			M5204L/Double	Cuff Width:14.5cm	70g	
03-0201-085	Adult Extra Large	White	M5104XL/Single	Cuff Length:50cm	80g	
03-0201-087			M5204XL/Double	Limb Circumference:25-35cm	80g	
03-0201-089	Large Adult	White	M5105/Single	Cuff Length:56cm	90g	
03-0201-091			M5205/Double	Limb Circumference:25-35/44cm	1pcs/Bag	
03-0201-093	Thigh	White	M5106/Single	Cuff Width:15.5cm	110g	
03-0201-095			M5206/Double	Cuff Length:69cm	120g	
				Limb Circumference:25-35/53cm	1pcs/Bag	
				Cuff Width:17.5cm	110g	
				Cuff Length:65cm	120g	
				Limb Circumference:34-47cm	1pcs/Bag	
				Cuff Width:22cm	160g	
				Cuff Length:78cm	170g	
				Limb Circumference:46-66cm	1pcs/Bag	

NIBP

### Reusable NIBP Cuff/Cone

(Polyester Cuff With TPU Bladder,Without Metal Loop)

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-149	Infant	White	M5141/Single	Cuff Width:6.5/4.0cm	25g	 <p>Tube Material:PVC Tube Length:20cm</p>
03-0201-151			M5241/Double	Cuff Length:21cm Limb Circumference:9-14cm	35g 1pcs/Bag	
03-0201-153	Pediatric	White	M5142/Single	Cuff Width:9.5/6.0cm	40g	
03-0201-155			M5242/Double	Cuff Length:32cm Limb Circumference:13-19cm	50g 1pcs/Bag	
03-0201-157	Small Adult	White	M5143/Single	Cuff Width:14/9.0cm	55g	
03-0201-159			M5243/Double	Cuff Length:44cm Limb Circumference:18-26cm	65g 1pcs/Bag	
03-0201-161	Adult	White	M5144/Single	Cuff Width:15/9.0cm	90g	
03-0201-163			M5244/Double	Cuff Length:50cm Limb Circumference:25-35cm	100g 1pcs/Bag	
03-0201-165	Large Adult	White	M5145/Single	Cuff Width:16.5/12cm	115g	
03-0201-167			M5245/Double	Cuff Length:64cm Limb Circumference:34-40cm	125g 1pcs/Bag	
03-0201-169	Thigh	White	M5146/Single	Cuff Width:21/14cm	170g	
03-0201-171			M5246/Double	Cuff Length:75.5cm Limb Circumference:39-50cm	180g 1pcs/Bag	

NIBP

### Cone Cuff With Metal Loop

(Polyester Cuff With TPU Bladder)

Used with: ABPM or Electronic BP Meter & Sphygmomanometer

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-129	Small Adult	White	M5133/Single	Cuff Width:12/9.0cm Cuff Length:39cm	90g 115g	 Tube Material:PVC Tube Length:50cm
03-0201-131			M5233/Double	Limb Circumference:18-26cm	1pcs/bag	
03-0201-133	Small Adult Extra		M5133L/Single	Cuff Width:12/9.0cm Cuff Length:45cm	100g 125g	
03-0201-135			M5233L/Double	Limb Circumference:18-26/33cm	1pcs/bag	
03-0201-137	Adult		M5134/Single	Cuff Width:14/10cm Cuff Length:52cm	130g 155g	
03-0201-139			M5234/Double	Limb Circumference:25-35cm	1pcs/bag	
03-0201-141	Adult Extra		M5134L/Single	Cuff Width:14/10cm Cuff Length:56cm	140g 165g	
03-0201-143			M5234L/Double	Limb Circumference:25-35/44cm	1pcs/bag	
03-0201-145	Large Adult	M5135/Single	Cuff Width:17/12cm Cuff Length:62cm	150g 175g		
03-0201-147		M5235/Double	Limb Circumference:35-44cm	1pcs/bag		


**NIBP**

### Ambulatory Blood Pressure Cuff/Fan

(PU Cuff With TPU Bladder,With Metal Loop)

Used with: Ambulatory Blood Pressure Meter-ABPM

## Fan Style

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-173	Small Adult	White	M5153/Single	Cuff Width:12cm	140g	 <p>Tube Material:PVC Tube Length:100cm</p>
03-0201-175			M5253/Double	Cuff Length:45cm	190g	
03-0201-177	Adult	White	M5154/Single	Limb Circumference:18-26cm	160g	
03-0201-179			M5254/Double	Cuff Width:14cm	210g	
03-0201-181	Large Adult	White	M5155/Single	Cuff Length:57.5cm	1pcs/bag	
03-0201-183			M5255/Double	Limb Circumference:25-35cm	205g	
				Cuff Width:17cm	255g	
				Cuff Length:72cm	1pcs/bag	
				Limb Circumference:34-47cm		

NIBP


## Disposable NIBP Cuff

P/N	Size	Printing Colour	Model	Description	Weight&Package	Photo
03-0201-185	Neonate-1#	Blue	M5541-1#/Single	Cuff Size:2.6x12.2cm	240g	 <p>Tube Material:PVC Tube Length:20cm 20pcs/Box</p>
03-0201-187			M5641-1#/Double	Limb Circumference:3-5.5cm	250g	
03-0201-189	Neonate-2#	Blue	M5541-2#/Single	Cuff Size:3.2x15.4cm	270g	
03-0201-191			M5641-2#/Double	Limb Circumference:4-8cm	280g	
03-0201-193	Neonate-3#	Blue	M5541-3#/Single	Cuff Size:4.2x18.7cm	310g	
03-0201-195			M5641-3#/Double	Limb Circumference:6-11cm	320g	
03-0201-197	Neonate-4#	Blue	M5541-4#/Single	Cuff Size:5.2x23cm	370g	
03-0201-199			M5641-4#/Double	Limb Circumference:7-13cm	380g	
03-0201-201	Neonate-5#	Blue	M5541-5#/Single	Cuff Size:6x25cm	410g	
03-0201-203			M5641-5#/Double	Limb Circumference:9-14.5cm	420g	
03-0201-205	Infant	Orange	M5542/Single	Cuff Size:5.7x21.5cm	240g	 <p>Tube Material:PVC Tube Length:20cm 1pcs/Bag</p>
03-0201-207			M5642/Double	Limb Circumference:9-14.5cm	280g	
03-0201-209	Pediatric	Green	M5543/Single	Cuff Size:8.2x30.5cm	26g	
03-0201-211			M5643/Double	Limb Circumference:13-21.5cm	36g	
03-0201-213	Small Adult	Sky Blue	M5544/Single	Cuff Size:10.5x41.5cm	37g	
03-0201-215			M5644/Double	Limb Circumference:21-27cm	47g	
03-0201-217	Adult	Sky Blue	M5545/Single	Cuff Size:13.5x51.5cm	60g	
03-0201-219			M5645/Double	Limb Circumference:26-35.5cm	70g	
03-0201-221	Adult Extra	Sky Blue	M5545L/Single	Cuff Size:13.5x66cm	60g	
03-0201-223			M5645L/Double	Limb Circumference:26-35.5/45cm	70g	
03-0201-229	Large Adult	Deeep Red	M5546/Single	Cuff Size:16.6x66.5cm	100g	
03-0201-231			M5646/Double	Limb Circumference:35-45cm	110g	
03-0201-233	Thigh	Gray	M5547/Single	Cuff Size:20.6x79cm	120g	
03-0201-235			M5647/Double	Limb Circumference:44-55cm	130g	

NIBP



### Animal Blood Pressure Cuff










P/N	Size	Printing Colour	Material	Description	Weight&Package	Photo
03-0201-237	Animal-1#	Blue	PVC Non-Woven	Cuff Size:2.6x12.2cm Limb Circumference:3-5.5cm	240g 20pcs/Box	 <p style="text-align: center;">Tube Material:PVC Tube Length:20cm</p>
03-0201-241	Animal-2#	Blue	PVC Non-Woven	Cuff Size:3.2x15.4cm Limb Circumference:4-8cm	270g 20pcs/Box	
03-0201-245	Animal-3#	Blue	PVC Non-Woven	Cuff Size:4.2x18.7cm Limb Circumference:6-11cm	310g 20pcs/Box	
03-0201-249	Animal-4#	Blue	PVC Non-Woven	Cuff Size:5.2x23cm Limb Circumference:7-13cm	370g 20pcs/Box	
03-0201-253	Animal-5#	Blue	PVC Non-Woven	Cuff Size:6x25cm Limb Circumference:9-14.5cm	240g 10pcs/Box	
03-0201-257	Animal-6#	White	Polyester	Cuff Width:8.8cm Cuff Length:29cm Limb Circumference:14-19cm	30g 1pcs/Bag	
03-0201-261	Animal-7#	White	Polyester	Cuff Width:11.5cm Cuff Length:39cm Limb Circumference:18-26cm	50g 1pcs/Bag	
03-0201-265	Animal-8#	White	Polyester	Cuff Width:14.5cm Cuff Length:50cm Limb Circumference:25-35cm	70g 1pcs/Bag	
03-0201-269	Animal-9#	White	Polyester	Cuff Width:17.5cm Cuff Length:65cm Limb Circumference:34-47cm	110g 1pcs/Bag	

NIBP

### NIBP Air Hose-TPU

P/N	Model	Use With Instrument	Weight&Package	Image
03-0202-001	M5146-147A Single	PHILIPS,Mindray 2.5m Long	110g 1pcs/Bag	
03-0202-042	M5146-173A Single	PHILIPS 2.5m Long	120g 1pcs/Bag	
03-0202-002	M5147-147A Single	PHILIPS-A1/A3,CSI Mindray,Spacelabs... 2.5m Long	110g 1pcs/Bag	
03-0202-003	M5148-147A Single	Create,Comen... 2.5m Long	110g 1pcs/Bag	
03-0202-004	M5149-151A Single	Criticare/CSI,Datascope Colin8800 Welch Allyn...	110g 1pcs/Bag	
03-0202-005	M5149-147A Single	2.5m Long	110g 1pcs/Bag	
03-0202-008	M5151-147A Single	Spacelabs(old)	110g 1pcs/Bag	
03-0202-009	M5154-153A Single	Colin 2.5cm length	110g 1pcs/Bag	






NIBP

P/N	Model	Use With Instrument	Weight&Package	Image
03-0202-011	M5157-147A Single	MEK 2.5m Long	110g 1pcs/Bag	
03-0202-012	M5169-147A Single	M&B,Pece Tech 2.5m Long	110g 1pcs/Bag	
03-0202-013	M5175-147A Single	Siemens/Drag	110g 1pcs/Bag	
03-0202-014	M5141-164A Double	GE-Marquette (Dash/Eagle) 2.5m long	210g 1pcs/Bag	
03-0202-041	M5141-172A Double	GE 2.5m long	220g 1pcs/Bag	
03-0202-046	M5142-164A Double	GE-Dinamap 2.5m long	210g 1pcs/Bag	
03-0202-016	M5142-157A Double	GE-Pro1000 2.5m long	210g 1pcs/Bag	
03-0202-017	M5143-163A Double	Nihon Konden 2.5m long	210g 1pcs/Bag	
03-0202-018	M5144-163/4A Double	GE-Datex 2.5m long	210g 1pcs/Bag	

NIBP

P/N	Model	Use With Instrument	Weight&Package	Image
03-0202-019	<b>M5146-155N</b> Neonate Single	PHILIPS,Mindray 2.5m long,Connect With:M5541-1#-4#	110g 1pcs/bag	
03-0202-020	<b>M5146-156N</b> Neonate Single	PHILIPS,Mindray 2.5m long,Connect With:M5541-1#-4#	110g 1pcs/bag	
03-0202-021	<b>M5147-156N</b> Neonate Single	PHILIPS,Mindray.... 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-022	<b>M5148-156N</b> Neonate Single	Create,Comen..... 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-023	<b>M5149-156N</b> Neonate Single	Spacelabs(old),CSl.. 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-026	<b>M5161-156N</b> Neonate Single	Spacelabs (new) 2.5m long,Connect With:M5541-1#-4#	110g 1pcs/bag	
03-0202-027	<b>M5154-156N</b> Neonate Single	Colin 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-047	<b>M5157-156N</b> Neonate Single	MEK 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-028	<b>M5169-156N</b> Neonate Single	M&B,Pace Tech 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	

**NIBP**










P/N	Model	Use With Instrument	Weight&Package	Image
03-0202-029	M5175-156N Neonate Single	Siemens/Drac 2.5m long,Connect with: M5541-1#-4#	110g 1pcs/bag	
03-0202-030	M5141-158N Neonate Double	Marquette/Dash 2.5m Long,Connect With:M5641-1#-4#	210g 1pcs/bag	
03-0202-031	M5142-158N Neonate Double	Pro1000 2.5m Long,Connect With:M5641-1#-4#	210g 1pcs/bag	
03-0202-045	M5143-158N Neonate Double	Nihon Konden 2.5m Long,Connect With:M5641-1#-4#	210g 1pcs/bag	
03-0202-033	M5144-158N Neonate Double	GE-Datex 2.5m Long,Connect With:M5641-1#-4#	210g 1pcs/bag	

**NIBP**










## Hose Adapter

P/N	Model	Discription	Weight	Image
03-0203-016	M163-147A	1-Hose Adapter 6cm Long	20g	
03-0203-001	M164-147A	1-Hose Adapter 6cm Long	20g	
03-0203-002	M165-147A	1-Hose Adapter 6cm Long	20g	
03-0203-017	M150-147A	1-Hose Adapter 6cm Long	20g	
03-0203-018	M152-147A	1-Hose Adapter 6cm Long	20g	
03-0203-006	M161-160-1A	1-Hose Adapter 6cm Long	15g	
03-0203-024	M163-160-1A	1-Hose Adapter 6cm Long	15g	
03-0203-007	M164-160-1A	1-Hose Adapter 6cm Long	15g	

NIBP







P/N	Model	Discription	Weight	Image
03-0203-008	M165-160-1A	1-Hose Adapter 6cm Long	15g	
03-0203-025	M150-160-1A	1-Hose Adapter 6cm Long	15g	
03-0203-026	M152-160-1A	1-Hose Adapter 6cm Long	15g	
03-0203-012	M163/163-147Y	1-Hose Y Adapter Connected with Single Cuff	25g	
03-0203-014	M163/164-147Y	1-Hose Y Adapter Connected with Single Cuff	25g	
03-0203-013	M165/165-147Y	1-Hose Y Adapter Connected with Single Cuff	25g	
03-0203-027	M163/163-160-1Y	1-Hose Y Adapter Connected with Single Cuff	20g	
03-0203-028	M163/164-160-1Y	1-Hose Y Adapter Connected with Single Cuff	20g	
03-0203-029	M165/165-160-1Y	1-Hose Y Adapter Connected with Single Cuff	20g	

NIBP

P/N	Model	Discription	Weight	Image
03-0203-030	M2T/2T-163Y	1-Hose Y Adapter Connected with Single Cuff	30g	
03-0203-023	M2T/2T-147Y	1-Hose Y Adapter Connected with Single Cuff	35g	
03-0203-031	M2T/2T-160-1Y	1-Hose Y Adapter Connected with Single Cuff	30g	
03-0203-015	M161-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	
03-0203-020	M163-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	
03-0203-032	M164-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	
03-0203-033	M165-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	
03-0203-034	M150-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	
03-0203-035	M152-147/147Y	2-Hose Y Adapter Connected with Double Cuff	35g	

NIBP



P/N	Model	Discription	Weight	Image
03-0203-036	M161-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	
03-0203-037	M163-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	
03-0203-038	M164-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	
03-0203-039	M165-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	
03-0203-040	M150-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	
03-0203-041	M152-157/157Y	2-Hose Y Adapter Connected with Double Cuff	30g	

NIBP

### Connector

P/N	02-0203-026	02-0203-049	02-0203-028	02-0203-029	02-0203-030	02-0203-050	02-0203-031
Model	CT-161	CT-161C	CT-162	CT-162N	CT-163	CT-163S	CT-163N
Photo							
Description	Cuff Connector Metal Male Bayonet PHILIPS Mindray...	Cuff Closure For closing CT-161 cuff connector during cleaning	Cuff Connector Plastic Male Bayonet PHILIPS Mindray...	Cuff Connector ID: 3/32" -2.4mm Male Bayonet For:Disposable Neonate Cuff	Cuff Connector ID: 5/32"-4.0mm Marqutte(2) Datex(1) Welch Allyn,Zoll...	Cuff Connector ID: 1/8"-3.2mm Marqutte(2) Datex(1) Welch Allyn,Zoll...	Cuff Connector ID: 3/32" -2.4mm For:Disposable Neonate Cuff
P/N	02-0203-032	02-0203-033	02-0203-034	02-0203-035	02-0203-055	02-0203-036	02-0203-037
Model	CT-164	CT-164S	CT-165	CT-165N	CT-165M	CT-166	CT-167
Photo							
Description	Cuff Connector ID: 5/32"-4.0mm Datex(1) Nihon Konden...	Cuff Connector ID: 1/8"-3.2mm Datex(1) Nihon Konden...	Cuff Connector Screw Pro1000(GE-Critikon) Welch Allyn...	Cuff Connector ID: 3/32" -2.4mm For: Disposable Neonate Cuff	Cuff Connector Metal Screw Pro1000(GE-Critikon) Welch Allyn...	Cuff Connector- For: PHILIPS Disposable Neonate Cuff	Cuff Connector Male Slip Luers For: Disposable Neonate Cuff

NIBP

P/N	02-0203-038	02-0203-056	02-0203-040	02-0203-041	02-0203-011	02-0203-012	02-0203-014
Model	CT-168	CT-168M	CT-170	CT-171-1	CT-150	CT-150N	CT-152
Photo							
Description	Cuff Connector Male Slip Luers For: Disposable Neonate Cuff	Cuff Connector Metal Male Slip Luers For: Disposable Neonate Cuff	Cuff Connector- PHILIPS, GE-Dinaclick	Cuff Connector- Universal	Cuff Connector- Female Luers CSI(old)...	Cuff Connector ID: 3/32" -2.4mm For: Disposable Neonate Cuff	Cuff Connector- Male Luers Spacelabs...
P/N	02-0203-015	02-0203-006	02-0203-007	02-0203-009	02-0203-010	02-0203-013	02-0203-017
Model	CT-152N	CT-146	CT-147	CT-148	CT-149	CT-151	CT-154
Photo							
Description	Cuff Connector ID: 3/32" -2.4mm For: Disposable Neonate Cuff	Monitor Plug PHILIPS Mindray Biolight...	Monitor Plug PHILIPS Mindray Spacelabs...	Monitor Plug Creative Comen...	Monitor Plug CSI(old)...	Monitor Plug Spacelabs...	Monitor Plug Colin

NIBP

P/N	02-0203-020	02-0203-021	02-0203-039	02-0203-043	02-0203-001	02-0203-002	02-0203-003
Model	CT-157	CT-157P	CT-169	CT-175	CT-141	CT-142	CT-143
Photo							
Description	Monitor Plug MEK Invivo...	Monitor Plug MEK Invivo...	Monitor Plug M&B...	Monitor Plug Siemens Drager	Monitor Plug GE (Dash,Solar...)	Monitor Plug GE (Pro1000...)	Monitor Plug Nihon Konden
P/N	02-0203-004	02-0203-016	02-0203-019	02-0203-022	02-0203-045	02-0203-051	02-0203-018
Model	CT-144	CT-153	CT-156	CT-158	CT-2T	CT-2TP	CT-155
Photo							
Description	Monitor Plug GE Dixtal	Connector >>CT-154 Colin	Connector >>CT-167, CT-168	Connector >>CT-167, CT-168	Connector	Connector	Connector >>CT-166 PHILIPS
P/N	02-0203-024	02-0203-044	02-0203-053	02-0203-054			
Model	CT-160-1	CT-Y	CT-172	CT-173			
Photo							
Description	Connector >>CT-171-1 Universal	Y Adaptor	Connector >>CT-170 GE	Connector >>CT-170 PHILIPS			

NIBP